# Septofil

FIL TRICOTÉ

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# Septofil

KNITTED CORD

Nonimpregnated tubular cord for gingival retraction

## Composition

Cotton cord.

## Indication

Retraction of the free gingiva before taking impressions.

## Directions for use

This cord produces a mechanical retraction. It makes a temporary opening of the sulcus possible, so as to take impressions or carry out some dental surgeries.

If gingival retraction is needed for several teeth, treat one tooth at a time.

Pull the cord out of the bottle and cut the needed length: the perimeter of the tooth level with the crown can vary from 2 to 6 or even 7 cm. Cut longer segments of cord than required so that when the cord is wound up around the tooth, the ends overlap. This makes the removal of the cord easier.

Place the cord in the sulcus (gingival crevice) using a metallic instrument with a blunt end so as not to injure the epithelial attachment. Start with the interproximal space because it is where the sulcus is usually deepest. Then, proceed clockwise until the complete tour of the tooth is done.

To obtain a greater retraction, once the cord has been inserted, it may be impregnated with a solution for retraction. Leave it in situ for the time indicated in the directions for use of the solution. The cord should be removed just before taking impressions except if the double impression technique is used. Slightly moisten the cord to make its removal easier and less painful.

## Warnings and precautions while using this medicine

- After impression-taking, the gingival sulcus should be inspected thoroughly to make sure that the epithelial attachment has not been torn and that no fragment of the retraction cord has been left in the gingival crevice, which would act as a foreign body and induce a severe local inflammation.
- This cord should only be used on healthy tissue, free from injury, bleeding and infection.

## Storage

Store the product in a dry place.

## Presentations

Bottle, with dispensing stopper, containing 2.5 meters of ultra-thin  
Bottle, with dispensing stopper, containing 2.5 meters of thin cord  
Bottle, with dispensing stopper, containing 2.5 meters of medium cord  
Bottle, with dispensing stopper, containing 2.5 meters of thick cord  
Box with 4 bottles, each of them containing 2.5 m of cord of assorted diameters (ultra-thin, thin, medium and thick).

*For professional use only In the practice of dentistry*

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